



FastGene®

Part compatibility between

FastGene® PAGE Protein System (PG01)

and

Bio-Rad Mini-PROTEAN Tetra, Handcast Gels 1.0 mm (#1658001EDU)

FastGene® PG01 → Bio-Rad #1658001EDU		
Does this FastGene® part fit into...	...this Bio-Rad part?	
PAGE components		
U sealing strip short for Biorad TGX gels (PG18)	Mini-PROTEAN Tetra Electrode Assembly (#1658037EDU) Mini-PROTEAN Tetra Companion Running Module (#1658038EDU)	YES
U sealing strip FastGene (PG19)	Mini-PROTEAN Tetra Electrode Assembly (#1658037EDU) Mini-PROTEAN Tetra Companion Running Module (#1658038EDU)	YES
U sealing strip long for 10 x 10 cm gels (PG20)	Mini-PROTEAN Tetra Electrode Assembly (#1658037EDU) Mini-PROTEAN Tetra Companion Running Module (#1658038EDU)	NO
Plastic dummy cassette short (PG23)	Mini-PROTEAN Tetra Electrode Assembly (#1658037EDU) Mini-PROTEAN Tetra Companion Running Module (#1658038EDU)	YES
Plastic dummy cassette long (PG13)	Mini-PROTEAN Tetra Electrode Assembly (#1658037EDU) Mini-PROTEAN Tetra Companion Running Module (#1658038EDU)	NO
Inner chamber with/without electrodes (PG11, PG12)	Buffer Tank (#1658039EDU)	NO
Gel/Blot chamber lid with electrodes	Buffer Tank (#1658039EDU)	YES
Gel casting components		
Glass plates short + long (1mm)	Gel casting frame	YES
Comb 1.0mm 10 wells/15 wells (part of PG10)	Glass plates short + long (1mm)	YES
Gel casting clip (PG22) with Glass plates short + long (1mm) inserted	Gel casting stand	NO
Sealing gaskets (PG07)	Gel casting stand	NO



Bio-Rad #1658001EDU → FastGene® PG01		
Does this Bio-Rad part fit into...	...this FastGene® part?	
PAGE components		
Sealing strip	Inner chamber with/without electrodes (PG11, PG12)	YES
Mini Cell Buffer Dam (#1653130EDU)	Inner chamber with/without electrodes (PG11, PG12)	YES
Mini-PROTEAN Tetra Electrode Assembly (#1658037EDU) Mini-PROTEAN Tetra Companion Running Module (#1658038EDU)	Gel/Blot chamber tank (PG05)	YES
Cell Lid With Power Cables (#1658041EDU)	Gel/Blot chamber tank (PG05)	NO
Gel casting components		
Glass plates short + long (1mm)	Gel casting clip (PG22)	YES
Comb 1.0mm 10 wells/15 wells	Glass plates short + long (1mm)	YES
Gel casting frame with Glass plates short + long (1mm) inserted	Gel casting base (PG21)	YES
Sealing gaskets	Gel casting base (PG21)	NO